



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: SILICA GEL NON-INDICATING

CAS No: 7631-86-9

Supplied as: Product code

CHE3194

Description

SILICA GEL NON-INDICATING 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Moisture at 145°C	%		2.0		
Water soluble sulphates	%		1.0		
Water soluble chlorides	%		0.050		

Material complies with BS 2540:1991 Granular Dessicant Silica Gel.

Revision date: 06/10/2020

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3259/MS/20201006