

Pharmacy Refrigerator

HYC-118A



LED Display

Microprocessor control, digital display and temperature adjustment with an increment of 0.1 °C.

Multiple Alarm Functions

Capable of alerting failures due to high and low temperature, sensor error, power failure, low battery, door ajar and remote alarm.

Ergonomic Design

810mm cabinet height makes it easy to place under test bench, saving storage space.

Lockable Design

Safety lock to prevent unauthorized access.

Superior Cabinet Temperature Uniformity, Auto Defrost

The test point of performance in the cabinet increased from 2 points to 9 points.

Energy Efficient Hydrocarbon (HC) Refrigeration System

Optimised refrigeration system, utilising high-efficiency HC compressor significantly reduces energy consumption.

Greener by Design

HC refrigerants are natural non-toxic and more environmentally friendly.

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Specifications



HYC-118A

Technical Data

| | |
|----------------|---------------------|
| Cabinet Type | Upright, glass door |
| Climate Class | N |
| Cooling Type | Forced air cooling |
| Defrost Mode | Auto |
| Refrigerant | HC |
| Noise (dB (A)) | 41 |

Performance

| | |
|-----------------|-----|
| Temp Range (°C) | 2-8 |
|-----------------|-----|



Control

| | |
|------------|----------------|
| Controller | Microprocessor |
| Display | LED |

Electrical Data

| | | |
|------------------------|---------------|-----------|
| Power Supply (V/Hz) | 220~240/50/60 | 115/50/60 |
| Power (W) | 230 | 230 |
| Electrical Current (A) | 1.5 | 1.5 |

Construction

| | | |
|----------------------------|---------|----------------|
| Capacity (L/Cu.Ft) | 118/4.2 | |
| Net/Gross Weight (approx) | (kg) | 46/51 |
| | (lbs) | 101.5/112.5 |
| Interior Dimension (W*D*H) | (mm) | 515*415*630 |
| | (in) | 20.3*16.3*24.8 |
| Exterior Dimension (W*D*H) | (mm) | 597*635*810 |
| | (in) | 23.5*25.0*31.9 |
| Packing Dimension (W*D*H) | (mm) | 680*690*910 |
| | (in) | 26.8*27.2*35.8 |

Accessories

| | |
|-----------------|-----|
| Caster | Y |
| Foot | Y |
| Porthole | Y |
| Shelves/Drawers | 3/1 |

Alarms

| | |
|---------------|---|
| High/Low Temp | Y |
| Remote Alarm | Y |
| Power Failure | Y |
| Sensor Error | Y |
| Low Battery | Y |
| Door Ajar | Y |

Certification

| | |
|----|---|
| CE | Y |
| UL | Y |