

SLS Lab Basics Pipettes				Pipette Tips				
SKU	Channel	Volume Range, μ L	Compatible Tips	Racked			Reload	Bulk
				Sterile, Filter	Sterile	Non Sterile		
SLS4900	Single	0.1 - 2.5	10 μ L, 10XL	SLE5146, SLE5148	SLE5112, SLE5110	SLE5124, SLE5122	SLE5136, SLE5134	SLE5234
SLS4902	Single	0.5 - 10	10 μ L, 10XL	SLE5146, SLE5148	SLE5112, SLE5110	SLE5124, SLE5122	SLE5136, SLE5134	SLE5234
SLS4910	Single	2 - 20	10XL	SLE5148	SLE5110	SLE5122	SLE5134	-
SLS4904	Single	5 - 50	20 μ L, 200 μ L*	SLE5150	SLE5114	SLE5126	SLE5138	SLE5236
SLS4906	Single	10 - 100	100 μ L, 200 μ L*	SLE5152	SLE5114	SLE5126	SLE5138	SLE5236
SLS4912	Single	20 - 200	200 μ L, 300 μ L	SLE5154, SLE5156	SLE5114, SLE5116	SLE5126, SLE5128	SLE5138, SLE5140	SLE5236
SLS4908	Single	100 - 1000	1000 μ L, 1250 μ L	SLE5158, SLE5160	SLE5118, SLE5120	SLE5130, SLE5132	SLE5142, SLE5144	SLE5236, SLE5242
SLS4914	Single	500 - 5000	-	-	-	-	-	SLS4042
SLS4916	Single	1000 - 10000	-	-	-	-	-	SLS4045
SLS4926	8	0.5 - 10	10 μ L, 10XL	SLE5146, SLE5148	SLE5112, SLE5110	SLE5124, SLE5122	SLE5136, SLE5134	SLE5234
SLS4918	8	5 - 50	20 μ L, 200 μ L*	SLE5150	SLE5114	SLE5126	SLE5138	SLE5236
SLS4920	8	10 - 100	100 μ L, 200 μ L*	SLE5152	SLE5114	SLE5126	SLE5138	SLE5236
SLS4922	8	20 - 200	200 μ L, 300 μ L	SLE5154, SLE5156	SLE5114, SLE5116	SLE5126, SLE5128	SLE5138, SLE5140	SLE5236
SLS4924	8	30 - 300	200 μ L, 300 μ L	SLE5154, SLE5156	SLE5114, SLE5116	SLE5126, SLE5128	SLE5138, SLE5140	SLE5236
SLS4936	12	0.5 - 10	10 μ L, 10XL	SLE5146, SLE5148	SLE5112, SLE5110	SLE5124, SLE5122	SLE5136, SLE5134	SLE5234
SLS4928	12	5 - 50	20 μ L, 200 μ L*	SLE5150	SLE5114	SLE5126	SLE5138	SLE5236
SLS4930	12	10 - 100	100 μ L, 200 μ L*	SLE5152	SLE5114	SLE5126	SLE5138	SLE5236
SLS4932	12	20 - 200	200 μ L, 300 μ L	SLE5154, SLE5156	SLE5114, SLE5116	SLE5126, SLE5128	SLE5138, SLE5140	SLE5236
SLS4934	12	30 - 300	200 μ L, 300 μ L	SLE5154, SLE5156	SLE5114, SLE5116	SLE5126, SLE5128	SLE5138, SLE5140	SLE5236

* Non-filter tips only