



Product Data Sheet

SLS8102/SLS8103

General Descriptions

Product Name	SLS Select centrifuge tube, 15mL
Further Description	Conical bottom, polypropylene (PP), white printed graduations with flat cap
General Applications	Manufactured from virgin ultra clear polypropylene to ensure good visibility Leak-resistant The cap grip assures easy opening and closing when wearing gloves White printed graduations with large marking strip for labelling
Container Material	Polypropylene
Cap Material	HDPE
Packaging	SLS8102 - Bagged SLS8103 - Recyclable card racks

General Features

Graduation Range	0.1 - 15.0mL
Graduation Interval	0.5mL
Temperature	-80 to 121°C
Bottom Graduation Range	0.1 - 1mL
Bottom Graduation Interval	0.1mL
DNase/RNase	Free
Non-Pyrogenic	<0.25EU/mL
Solvent Resistant Printing	Will not rub off
Leak-Proof at	-95kPa (14psi)
Max RCF	15,000xg <small>Refers to liquids with a specific gravity of 1.0, using appropriate rotors, cushions, and adapters</small>

