

26s gauge, Small Hub RN Needle, 2 in, point style 3, 6/PK

Specification Sheet (Part/REF # 7768-01)



Small Hub Removable (RN) needles can be used with removable (RN)-style syringes that are 100 μ L and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	26s gauge
Point Style	3
Needle Length	2 inch (51 mm)
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/needles/7768-01>