

22s gauge, Metal Hub Needle, 2 in, point style 2, 6/PK

Specification Sheet (Part/REF # 90038)



Metal hub (nickel plated brass) needles can be used with TLL syringes and LT or TLL connectors. The needles are 304 stainless steel. Standard needle length is 2" / 51 mm. Metal hub (N) needles are available in convenient six packs. Custom needle lengths and point styles are available.

Product Specifications

Gauge	22s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Metal (N)

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/needles/90038>