

Lovibond® Water Testing Tintometer® Group



LC100 SV100 Kit

Hand-Held Imaging Spectrocolorimeter With Viewing Stand



- Unique Housing Stand to Protect Sample from Ambient Light
- Enables Measurement of Food, Powders, Cosmetics & Pastes
- Versatile Communication of Colour Measurement
- Single, Simple, Reliable, Integrated Tool - No Need to Connect to PC

Part Number: 401120

Portable, Ergonomic & Extremely Versatile

The LC 100 / SV 100 combination gives you great flexibility in colour measurement. It is a lightweight tool that is easy to handle and ready to use. In just over a second, it not only provides you a reading of the sample but also a report telling how close the colour is to a measured standard. It is the ideal solution for the colour measurement of food, powder, cosmetics, pastes and pharmaceuticals.

Easy, Instant Reporting to Share and Store Results

The LC 100 / SV 100 gives you immediate results on the screen with comparisons to your set of standards. If you need a more formal, in-depth report, the LC 100 creates PDF and CSV files which can be accessed via a PC with the USB mass storage mode.

Enables Measurement of Food, Powders, Cosmetics & Pastes

The Lovibond® LC 100 / SV 100 integrated package offers you a unique opportunity to measure awkward substances such as food, powders, cosmetics and pastes. The SV 100 holder supports the LC 100 at an ergonomic viewing angle to enable the measurements in Lovibond® optical cells. A secure cap is provided to cover the sample, protecting it from any excessive ambient light, enabling quick yet accurate readings.

Versatile Communication of Colour Measurement

You can get results in just over a second and view them easily with:
Graphical Pass/Fail indicator
Display of ΔE for all common methods ie CIELAB, CMC, CIE 94, 2000
Display of measured and delta differences with $L^*a^*b^*C^*h^\circ$ for standard and sample
Graphical $L^*a^*b^*$ plot
Verbal colour difference description

Single, Simple, Reliable, Integrated Tool - No Need to Connect to PC

Utilising the LC 100 / SV 100, you can ensure that all your staff are speaking the same language with an easy pass/fail measurement.

Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Oil and Gas | Pharmaceutical Industry

Application

Chemicals | Coatings | Food and Beverage | Household Products | Medical and Clinical Products | Pharmaceuticals and Cosmetics | Plastics | Sugar Solutions, Syrups and Honey | Textiles

LC100 SV100 Kit

A very easy to use, extremely portable instrument spectrophotometer for quick yet accurate reporting of any possible deviations of colour against your defined colour scale. Up to 350 sample measurements can be stored on the device and compared to reference standards.

The SV 100 holder provide even more flexibility enabling you to measure samples of food, powders, cosmetics and pastes.

Delivery Scope Text

- Sturdy Case Including
- LC 100 instrument
- Carrying pouch
- 10 x 10 mm plastic cells
- USB cable for power supply and data transfer
- Viewing stand
- Cell adaptor

Technical Data

Applicable Standard	ASTM E308
Colour Type	Transparent Translucent Opaque
Interfaces	USB
Illuminants	CIE Illuminant A CIE Illuminant D65
Portability	Portable

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
Fax: +49 (0)231/94510-30
sales@tintometer.de
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
water.sales@tintometer.com
www.lovibond.com
UK

Tintometer AG

Hauptstraße 2
5212 Hausen AG
Tel.: +41 (0)56/4422829
Fax: +41 (0)56/4424121
info@tintometer.ch
www.tintometer.ch
Switzerland

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebu Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@tintometer.com.br
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@tintometer.us
www.lovibond.com
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@tintometer.com
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain